# THE "RECOMBINANT PROTEINS" TECHNICAL PLATFORM

## IS A TECHNICAL DEPARTMENT OF INEM DEDICATED TO THE PURIFICATION OF RECOMBINANT PROTEINS BY CHROMATOGRAPHY

This platform provides some equipment to allow you to carry out the production and purification of your recombinant protein, or to fractionate biological samples.

The "Recombinant proteins" technical platform is in a cold room located on **the 1st floor** of the **Necker Faculty of Medicine building**, which houses the Institut Necker Enfants Malades (INEM).

### LOCATION AND CONTACT

This technical platform is under the responsibility of:

- A part-time operational manager: Emeline Pacreau, assistant engineer  
Contact: [emeline.pacreau@inserm.fr](mailto:emeline.pacreau@inserm.fr) Tel: 01.40.61.55.06
- A Scientific Manager: Pamela Schnupf, INSERM research director  
Contact: [pamela.schnupf@inserm.fr](mailto:pamela.schnupf@inserm.fr)


### ACCESSIBLE EQUIPMENT

The material listed below is made available to users; the maintenance of the devices is for the moment ensured by the resources INEM. The purchase of consumables that are not available are the responsibility of the users (eg: disposable columns, buffers, dialysis cassettes, etc.).

- An **AKTA-FPLC PURIFICATION SYSTEM** (chromatography system with high pressure)

This system is controlled by Unicorn Software version 5.3 installed on a laptop.

A few disposable columns (ion exchange, "His6" type affinity) can be made available to users according to availability.

- Alongside the purification of recombinant proteins, chromatography equipment can also be used to fractionate biological samples (serum, tissue lysate, conditioned culture media, etc.).

- A **FRENCH PRESS**, a device that lyses bacteria at high pressure (~1,000 psi).

This step releases a soluble fraction and an insoluble fraction. The recombinant proteins over-expressed in *E. coli* are very often found in the insoluble part (inclusion body), which facilitates the subsequent purification steps but requires a suitable solubilization/renaturation protocol (i.e. "refolding").

This equipment is suitable for the lysis of large volume cultures (example: 1L of culture resuspended in 40 mL of buffer). It is not suitable for lysing small volumes of cultures.

- A **CONTINUOUS DIALYSIS SYSTEM** for the renaturation of proteins produced in insoluble form (inclusion bodies).

Insoluble proteins must first be dissolved in denaturing buffers (8M urea, 6M guanidium chloride) then renatured by progressive elimination of these denaturing agents during continuous dialysis. The "home-made" system includes a 5-liter tank equipped with an overflow, in which are placed the dialysis cassettes loaded with the protein solution to be renatured, an 8- liter fresh buffer tank, and a peristaltic pump allowing continuous dialysis.

## PROCEDURE FOR USE

Anyone interested in using the platform will contact the managers in order to present their scientific project. Members of the platform are available to discuss upstream protocols and downstream results (see Discussion sheet).

The production steps upstream of the purification (bacterial cultures, verification of the production by electrophoresis, etc.) will be carried out by the users within their team (bacteriology room of the INEM available). The manipulation of the purification system will be ensured by the Assistant Engineer of the Technical Platform, except for recurring users who will follow an ad hoc training if they wish to become autonomous. Access to the platform will remain controlled.

## RECOMBINANT PROTEIN PLATFORM

### DISCUSSION SHEET

1- Has your protein of interest already been produced (published)?

- a. Yes: we look together, and we follow the same protocol
- b. No: see point 2

2- Does the expression plasmid exist?

- a. Yes: we attempt a "classic" production test – see point 3
- b. No:
  - i. Reflect together on a suitable expression plasmid (tags, etc.)
  - ii. You (not the platform) builds, sequences, and amplifies the vector

Nicolas Kuperwasser (SFR Genome editing platform) can probably help you

### 3- TESTS of analytical production:

- a. We discuss the protocol together, we help you if necessary
- b. First test in small volume (10-50 ml)
- c. SDS PAGE gel before/after induction, reducing/non-reducing (if relevant)
- d. Production analysis (to be done together)
  - i. If OK: see point 4
  - ii. If no production: we think ...

### 4- Production of a first preparatory batch

- a. 500-1000 ml – standard protocol or adapted according to your protein/protocol
  - i. At home (agitator available bacterial room) – we help you if needed
  - ii. Check production by SDS PAGE (same as analytical)
- b. IF OK, French press (breaking of bacteria)
  - i. By the platform (autonomy possible if recurring use)
  - ii. Separate soluble/insoluble fraction
  - iii. By you: SDS PAGE to determine if soluble or insoluble protein
    1. If insoluble: see point 5
    2. If soluble: see point 6

### 5- If insoluble protein:

- a. Standard denaturation protocol (urea/continuous dialysis) – we show you the first time
  - i. If protein precipitates during renaturation: we reflect...
  - ii. If protein remains soluble: proceed to purification by chromatography

### 6- Purification by chromatography

- a. Depending on the characteristics of what needs to be purified (volume, pH, pI of the protein, etc.), we decide on a type of column (Emeline does the purification, at least at the beginning):
  - i. Ion exchange (Q or S): above (Q) or below (S) the pI of the protein of interest
  - ii. Affinity: Tag His 6, ...
  - iii. Sieve (gel filtration): analytical test on AKTA, preparative on Gradifrac
- b. For each purification, you make an SDS PAGE of the fractions/peaks to locate the protein, define its state of purity etc.
  - i. If OK, you keep it at home (frozen, freeze-dried), then proceed to the functional tests
  - ii. If not OK: think about purification protocol, change column, etc.

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